

BABY BOARD

Pediatric immobilization and trasport

Code: 0313922









www.ferno.it

BABY BOARD





The **Medkids Baby Board** model 676 helps maintain cervical neutral alignment in infants with the self-contained **Pneumatic Positioning Device** (*PPD*). It is used to restrain, position and immobilize an infant during the transport in an incubator or on an ambulance cot.

MAIN FEATURES

- Made of durable Dupont Cordura and it is X-ray translucent.
- Compatible with most EMS cots and isolettes to provide safe and secure transport.
- Crescent-shaped head blocks secure cervical spines for infants too small for standard size cervical collars.
- Unique design accommodates all central lines, umbilical lines, chest tubes, urinary catheters and other invasive and non-invasive equipment.
- Designed for infants ranging from 2 lb. (1 kg) to 15 lb. (7 kg).
- Can be placed in a transport incubator.
- Simple care and maintenance.
- Easily-replaceable parts.

SPECIFICATIONS

Length	580 mm
Width	280 mm
Height	30 mm
Weight	2,5 kg
Load capacity	2-15 lb (1 - 7 kg)



