

# TSL Expander Kit

*Expander for Scoop EXL*

Codice: 25-0030-002



# TSL Expander Kit



**Scoop EXL code: 0708002**

**TSL Expander** is an accessory made from high-strength aluminum alloy that attaches to the Scoop EXL atraumatic stretcher and its standard TSL hook fitting.

Designed for the movement and transport of a large patient in all intra- and extra-hospital environments.

The option of using the Expander Kit, the Scoop EXL atraumatic stretcher, and Easy Lift hangers, also allows you to use the system with lifters.

The kit contains: 2 x Expanders (red), 3 x restraint extensions, and 1 x container bag.

## KEY SPECIFICATIONS

- Kit contents can be stowed anywhere.
- High-strength metal restraints
- Assembles in seconds.
- Fully decontaminable.
- Mortuary belts available on request
- Compatible with Easy Lift hangers



**Code: 0313659-IT-M**



**Code: PS0765EXL**

## OPTIONAL ACCESORIES

Description	Code
TSL Expander Kit	25-0030-002
Easy Lift (optional)	PS0765EXL
Mortuary belts (optional)	0313659-IT-M



## TECHNICAL SPECIFICATIONS

Height	150 mm
Width	300 mm
Weight (KIT)	1.5 kg
<b>Load capacity</b>	<b>227kg</b>