## HARMONIZED AND NON-HARMONIZED STANDARDS



- EN 1498:2006 Human Lift
- EASA CM No.: CM-CS-005 issue 01 Appendix 1, CS 27/29.865 Human external cargo (HEC)
- EN ISO 14971:2019/A11:2021 Risk assessment
- EN ISO 15223-1:2021 Labeling
- EN ISO 13485:2016+A11:2021 Medical quality
- EN ISO 9001:2015 Quality
- CE Marking
- Complies with the essential requirements listed in Annex I of the Regulation (EU) 2017/745 concerning Medical Devices

cod. XT PRO MILITARY
NATO number 6515150288871



## **EQUIPMENT**

- n°1 Board made of composite material
- n°1 Carry bag
- n°1 Qhi
- n°1 Triangular head immobiliser
- n°1 Chin strap
- n°4 Restraints
- n°2 Auxiliary handles
- n°2 Padded leg hoops
- n°1 Snap hook
- n°1 Lifting tape
- n°1 Filling cushion

## **GFERNO**°

Italy
P +39 051 6860028
F +39 051 6861508
info.it@ferno.com

www.ferno.it

Schweiz
P +41 41 259 60 00
info.ch@ferno.com
www.ferno-schweiz.ch



