

NOTICE OF CHANGE

PacRac+ Equipment Table

Doc ID: 00435_NOC_PacRac+ IT



The PacRac+® Equipment Table is an accessory that mounts to an ambulance cot or iN|X®. The PacRac+ is designed to secure medical instruments such as a defibrillator and/or heart monitor.

Effective from September 2018, the Model FWEPR+ PacRac+ Equipment Table will be upgraded, and an additional model for use with the POWER X1 added to the PacRac+ range.

The former FWEPR+ table has been superseded with a new table top featuring an integrated Ferno iN|TRAXX system within the table's surface. All other aspects of the PacRac+ remain unchanged, including the provision of equipment restraints.

As the global adoption of the iN|TRAXX Integrated Vehicle Component System continues to grow, we are confident the upgraded table top will future-proof and extend the versatility of the PacRac+ as well as provide users with an additional fastening option during patient transport.

The PacRac+ is compatible with Ferno® Mondial® Monobloc series cots, POWERX (UK), POWER X1 (US) and with the iN|X® Integrated Patient Transport & Loading System™. Model FWEPR+ IT for POWER X1 (US) is 1 in / 25 mm wider than other models to accommodate the cot's wider width. Detailed specifications are available from Ferno Australia's **website**.

PacRac+ Model	Compatible with	Additional Mount Kits required
FWEPR+ IT	<ul style="list-style-type: none"> iN X Mondial Monobloc 	Included with PacRac+
FWEPR+ IT	<ul style="list-style-type: none"> POWERX (UK) 	Fixing Kit PacRac+ for POWERX (009010150)
FWEPR+ PX1	<ul style="list-style-type: none"> POWER X1 (US) 	Included with PacRac+