

# WIZLOC JUNGLE

*Multi-sized Pediatric Cervical Collar*

Code Number 21-0140-020



Rel. 110116

# WIZLOC JUNGLE



Ferno's WizLoc JUNGLE collar is a multisized cervical immobilizer specifically studied for traumatized patients in the age bracket including 2 to 10 years. The pediatric Wizloc Jungle collar allows you to manage all types of injuries resulting from head trauma.

## FEATURES

- Mono-valve use: easy and quick positioning.
- Bi-valve use: can be adapted to patients with particular needs.
- Can be placed into position either from the left or right side of the patient according to the patient or operator needs.
- Micrometric front adjustment: major adaptability to the patient and neutral spine alignment.
- Rear section adjustment: can be used on an ample range of pediatric patients without occipital compression.
- Adjustment directly on the patient with collar in place.
- Rear section can easily be shaped and adapted to any form of the occipital area.
- Hypoallergenic lining, closed cell, to avoid absorption of biological liquids.
- Totally x-ray transparent and MRI compatible.
- Can be totally decontaminated (up to 10 times).
- Numbered decontamination stickers to affix directly on the collar (provided).
- Stickers and album (standard equipment) to entertain and comfort the patient.
- Behavior instructions for parents translated into the most common languages spoken.
- Ample range of use: pediatric patients from 2 to 10 years of age.
- Training course on pediatric trauma management from the Ferno Academy.

## TECHNICAL SPECIFICATIONS

Max Length	open: 450 mm
Height: front	160 mm
Height: back	150 mm
Weight	180 g
Range of application	2-10 years old

### Code

**21-0140-020**

### Description

WizLoc JUNGLE, multisized pediatric cervical collar, comes with adhesive numbers to keep track of how many times it has been decontaminated and with an album and stickers for children.