


XT LINE

Extrication devices



XT Plus

XT part no. Kit XT PLUS-B

CONTENT: no. 1 Carbon board, no. 1 Set of classic restraints, no. 1 QHI, no. 1 Carry bag,
no. 1 Set of Wiz Loc collars (3 Pcs), no. 1 Triangular head immobiliser



XT Plus

Ferno XT: Very lightweight, easy to carry and user-friendly for any kind of operator, it allows a complete, easy and quick immobilization of the cervical spine.

XT Plus part no. Kit XT PLUS-B



no. 1 carbon board
XT CARBON BOARD
part no. 10-1600-001

no. 1 Set of restraints
SET OF CLASSIC RESTRAINTS
part no. Yellow buckle. 10-1602-002
part no. Black buckle. 10-1604-002
part no. Green buckle. 10-1603-002
part no. Red buckle. 10-1601-002

no.1 QHI, triangular head immobiliser,
neoprene strap
**QHI, TRIANGULAR HEAD IMMOBILISER,
NEOPRENE STRAP**
part no. QHI-B-N



no. 1 carry bag
XT CARRY BAG
part no. 10-1605-001

no. 1 Set of Wiz Loc collars
CONTAINING 3 PCS
part no. 0822074

MAIN FEATURES

- Designed for spine immobilization and for the **safe extrication of patients from vehicles or small spaces**.
- It can be used for technical rescue operations in difficult environments (water and mountains) thanks to the Pro and Floating kits.
- Designed to perform **quick extrication** and **precautionary extrication** operations in a safe way for the operator.
- The **4 safety restraints**, two on the chest in crossed configuration and two on the groin, allow to immobilise the patient without compressions, ensuring the maximum comfort and safety.
The buckles have different colours for immediate identification and fast application.
- The **4 red lifting and carrying handles** allow a steady and ergonomic grip by the operator and ensure an efficient extrication of the patient.
- The **QHI-B Quick Head Immobilizer** guarantees an effective head-cervical spine immobilisation, thanks to the combination with the neoprene triangular head immobiliser.
- Extremely lightweight thanks to its structure in carbon composite fibre.
- Always ready for use, it comes in a convenient carry bag.
- Easy to sanitise.
- Compatible with non-traumatic stretchers, spinal boards and vacuum mattresses.

NATO NUMBER:

XT COMPLETE-B
6530150200590

XT PRO
6530150200591

XT FLOATING
6530150200592

XT

WARNING: If you already have an XT, you can update it to the PRO version thanks to the special kit.

XT Pro part no. XT PRO

CONTENT: no. 1 carbon board, no. 1 Kit Pro, no. 1 QHI, no. 1 carry bag, no. 1 Triangular head immobiliser.

The human lift kit allows vertical rescue in small spaces.

The device keeps the patient posture on the rescue line, considerably reducing the handling space.

Configuration speed is ensured by the single tape that facilitates patient lifting.

The Quick Head Immobilizer (QHI-B) allows the operator to efficiently immobilise the head-cervical spine area of the patient



XT part no. XT PRO



no. 1 carbon board
XT CARBON BOARD
part no. 10-1600-001



no. 1 Kit Pro
CONTENT:
no. 4 restraints with Cobra buckle,
no. 2 handling rings,
no. 2 paddings for restraints
(leg hoops),
no. 1 snap hook,
no. 1 lifting tape
part no. KIT PRO



no. 1 carry bag
XT CARRY BAG
part no. 10-1605-001



no.1 QHI, triangular head
immobiliser, neoprene strap
**QHI, TRIANGULAR HEAD
IMMOBILISER, NEOPRENE STRAP**
part no. QHI-B-N

Human LIFT



- Restraint system tested in accordance with EN1498
- Restraints with Cobra buckles - capacity 9 kN (917 kgf) made of aluminium alloy 7075
- Hanging ring - capacity 45 kN
- Leg hoops to ensure greater comfort

XT

XT Floating part no. XT FLOATING

CONTENT: no. 1 carbon board, no. 1 Kit floating, no. 1 QHI, no. 1 Triangular head immobiliser, no. 1 Special carry bag.

Xt floating is a kit designed to ensure extreme manoeuvrability and reliability during water rescue.

The head float and the two side floats ensure an optimal patient support, preventing capsizing and protecting him or her from accidental impacts.

The layout of semi-rigid polyurethane floating tubes is conceived to **keep the patient with the head out of water.**

Buoyancy is greater than the one of a buoyancy vest. The **FAST CLIP** technology allows to fasten the floating system in a few seconds, thus considerably speeding up water rescue operations.

The Quick Head Immobilizer (QHI-B) allows the operator to efficiently immobilise the head-cervical spine area of the patient.



XT part no. XT FLOATING



no. 1 carbon board
XT CARBON BOARD
part no. 10-1600-001



no. 1 Kit Floating
CONTENT:
no. 4 handling rings,
no. 3 floating tubes,
no. 2 transverse retaining strips
(Yellow and Red),
no. 1 Longitudinal Retaining Strip


part no. KIT FLOATING



no. 1 carry bag
**XT CARRY BAG
FLOATING**
part no. 10-1605-002



no.1 QHI, triangular head
immobiliser, neoprene strap
**QHI, TRIANGULAR HEAD
IMMOBILISER, NEOPRENE STRAP**
part no. QHI-B



XT

PART NO.

DESCRIPTION

XT PLUS-B

XT
no. 1 carbon board, no. 1 Set of restraints classiche, no. 1 QHI,
no. 1 Triangular head immobiliser, no. 1 carry bag,
no. 1 Set of Wizloc collars (3 pcs)

XT PRO

XT + KIT PRO
no. 1 carbon board, no. 1 Set of restraints pro, no. 1 QHI,
no. 1 carry bag, no. 1 Triangular head immobiliser

XT FLOATING

XT + KIT FLOATING
no. 1 carbon board, no. 1 Set of restraints floating, no. 1 QHI,
no. 1 Triangular head immobiliser, no. 1 carry bag floating

WARNING: If you already have an XT, you can update it to the PRO or FLOATING version thanks to the special kits.

KIT PRO




KIT PRO

no. 4 restraints with Cobra buckle, no. 2 handling rings,
no. 2 paddings for restraints (leg hoops), no. 1 snap hook,
no. 1 lifting tape

PART NO. KIT PRO

KIT FLOATING



KIT FLOATING

no. 4 handling rings, no. 3 floating tubes,
no. 2 transverse retaining strips (Yellow and Red),
no. 1 longitudinal retaining strip

PART NO. KIT FLOATING

XT

SPARE PARTS

Description	Part no.
Head immobiliser with triangle and neck strap (black bag included)	QHI-B-N
Chin strap for QHI	25-0601-002
Neoprene triangular head immobiliser	038245211
Set of 4 restraints, plastic buckles (XT standard)	25-0500-00
Storage and carry bag (XT & XT Pro)	10-1605-001
Storage and carry bag (Floating)	10-1605-002

TECHNICAL SPECIFICATIONS XT

Length	830 mm
Width	300 mm
Height	60 mm
Weight	2 kg
Load capacity	160 kg

TECHNICAL SPECIFICATIONS - XT Floating (XT + Floating kit)

Length	1200 mm
Width	350 mm
Height	110 mm
Weight	3,3 kg
Load capacity	160 kg

TECHNICAL SPECIFICATIONS - XT Pro (XT + PRO kit)

Length	830 mm
Width	300 mm
Height	60 mm
Weight	3,4 kg
Load capacity	160 kg



Italy
T: +39 051 68 600 28
info.it@ferno.com

Switzerland
T: +41 41 259 60 00
info.ch@ferno.com

